Agenda

- Start Measure a Bean lab
- Class time to work on study guide



Measure a bean objective

- The purpose of this lab is for you to become familiar with the <u>names</u> and <u>uses</u> of the laboratory equipment.
- Record all measurements and data on the <u>data table</u> on <u>page 4</u>.




Day one materials per lab group

- 1. Five beans in a paper cup
- 2. One sharpie marker
- 3. One classroom ruler
- 4. One beaker of water
- 5. One graduated cylinder
- 6. One paper towel
- 7. One electronic balance (shared among the entire class)



Day one procedure

- 1. Number 5 beans 1-5 with a sharpie.
- 2. Measure the length and width of *each bean in mm* with a metric ruler.
- 3. Use beaker to fill a graduated cylinder to the 10 mL mark
- 4. Add *all 5 beans* to the graduated cylinder. To determine the volume of the beans, see how many mL the water raised. Don't forget to subtract the 10 mL you started with!
- 5. Dry your beans with a paper towel. Record the mass of *all 5 beans* using an electronic balance.
- 6. Write your team name on the cup with a sharpie, place *all five beans* in the cup, and fil it with 50 mL of water.
- 7. Place your cup of beans and water on the back table.
- 8. Return your lab equipment where you got it and clean up your area.
- 9. Work on the day one lab questions.



Data tables

DATA TABLES (20 points, 10 points each)


	DAY ONE		DAY TWO	
Bean #	Width (mm)	Length (mm)	Width (mm)	Length (mm)
1				
2		2		
3				
4				~
5				s.
Mean				_
	1			
Volume of				
all 5 beans				
(mL)				
Mass of all	0		2	
5 beans (g)				
10	2		<i>8</i>	

How to use the electronic balance

- 1. Make sure the balance is on
- 2. Press the On/Zero button and wait until you see 0.0 on the screen
- 3. Place your beans on the balance
- 4. Be careful not to lean on the table
- 5. Record your data



READ THE DAY ONE PROCEDURE!!!!!



FINISHED WITH THE DAY 1 QUESTIONS AND CLEANED UP???

Work on your study guide
Metric system word searches
Quiet white board review